



	No	Type
AB Dental		
Narrow Platform (NP)	55021	11
Internal Hex Platform (SP)	55049	41
Conical Platform (CP)	55046	38
Multi-Unit Abutment	55080	72
ADIN		
S/OS Touareg	55034	26
UNP Touareg Closefit 2.75	55059	51
NP Touareg Closefit 3.0	55059	51
RP Touareg CloseFit 3.5	55046	38
WP Touareg CloseFit 4.3, 5.0	55020	10
Swell 3.3	55034	26
UniFit 3.5-6.0	55066	58
AIRES		
3.3	55049	41
Alfa Gate		
Bioactive All	26	55034
MAX All	26	55034
NP M+ 3.5	36	55044
RP M+ 4.3-5.0	29	55037
Porous All	26	55034
S-Line 3.0	51	55059
Alpha Dent		
Narrow Platform Conical connection	55025	16
Wide Platform Conical connection	55057	49
Narrow Platform Hexagon connection	55039	31
Regular Platform Hexagon connection	55057	49
Alpha-Bio		
Internal Hex	55034	26
Conical Hex 3.2	55065	57
Conical Hex 3.75-5.0	55051	43
Andrew Medical		
Conical Type K Lasermade All diameters	55066	58
Anker		
SB III 3.5	55023	14
SB III 4.0-5.0	55020	10

	No	Type
Ankylos		
C/X	55025	16
Anthogyr		
Axiom Bone Level 3.4 - 5.2	55065	57
Axiom Trilobe	55065	57
Narrow, Regular Axiom Tissue Level	55050	42
Argon		
Rapid 2mm K3pro 3.0 - 3.5	55026	17
Rapid 3mm K3pro 4.0 - 6.0	55019	9
Rapid short 3mm K3pro 4.0 - 6.0	55044	26
Astra Tech		
EV Osseospeed 3.0	55058	50
EV Osseospeed 3.6	55046	38
EV Osseospeed 4.2	55057	49
EV Osseospeed 4.8	55055	47
EV Osseospeed 5.4	55048	40
EV Profile Osseospeed 4.2	55034	26
TX Osseospeed 3.0 S	55064	56
TX Osseospeed 3.5-4.0 S	55046	38
TX Osseospeed 4.5-5.0 S	55017	7
Avinent		
External Coral & Ocean 3.5	55038	30
External Coral & Ocean 4.1	55011	1
External Coral 5.1	55013	3
Internal Coral & Ocean 3.5, 4.1	55051	43
Conical Ocean 3.5	55066	58
Conical Ocean 4.1	55004	18
B&B Implant		
Duravit Slim	55021	11
Duravit 3P	55035	27
Duravit Evolution	55035	27

	No	Type
Dentegris		
3.3	55034	26
3.75	55057	49
4.5	55049	41
3.8/4.1	55049	41
5.5	55020	10
Dentis		
Cleanlant S-Clean 3.7-4.8	55020	10
Narrow S-Clean OneQ 3.3	55064	56
Dentium		
Superline	55037	29
NR Line	55066	58
Simple Line II GH 1, 6	55014	4
Dentsply		
PrimeTaper/OmniTaper 3.0	55058	50
PrimeTaper/OmniTaper 3.6	55046	38
PrimeTaper/OmniTaper 4.2	55057	49
PrimeTaper/OmniTaper 4.8	55055	47
PrimeTaper/OmniTaper 5.4	55055	47
DEVA		
Bone Level 3.5	55064	56
Bone Level 4.5	55051	43
DIO		
Submerged SM 3.8	55049	41
Submerged SM 5.0	55020	10
Narrow UF II	55064	56
Regular UF II	55020	10
Ditron Dental		
MPI 3.3-6.0	55087	79
Ulitimate 3.75-6.0	55087	79
DSI		
Regular Internal Hex	55049	41
Narrow Internal Hex	55039	31
Regular Conical connection	55037	29
Narrow Conical connection	55046	38

	No	Type
DSP Biomedical		
External Standard (HE) Hexagonal	55014	4
Internal Standard (HiS) Hexagonal	55034	26
Internal Morse (HIM) Hexagonal	55046	38
Internal Large (HiL) Hexagonal	55037	29
Morse Taper (CMI) Indexed	55025	16
DYNA Dental		
Helix DC	26	55034
Helix ST	26	55034
Helix TM	16	55025
Easy Implant		
Ocean	55046	38
Lilas	55037	29
Edierre		
3-lobe	55021	11
4-lobe	55025	16
Edison Medical		
MUA	55033	25
Euroteknika		
Naturactis 3.5-5.0	55081	73
Naturactis/Naturall+ 3.0	55064	56
Uneva	55011	1
Futura implants		
PRE 4.0	55014	4
PRI 4.0	55037	29
galimplant		
ICI 3.2 & 3.8	55046	38
iPx 4012	55046	38
MUSA	55128	94
GC Implant		
Narrow Aadva	55081	73
Regular Aadva	55066	58
Wide Aadva	55066	58
Short Aadva 4.2	55044	36
Plus Genesio 3.4-5.0	55085	77
Plus SETiO 3.8	55011	1

	No	Type
Keystone		
I-HEXMRT 7.0-9.0	55051	43
Narrow Prima/Connex	55049	41
Regular Prima/Connex	55067	59
Wide Prima/Connex	55041	33
Narrow Prima/Genesis	55049	41
Regular Prima/Genesis	55067	59
Wide Prima/Genesis	55041	33
Narrow Restore	55079	71
Regular Restore	55011	1
Wide Restore	55013	3
TILOBEMAXX 7.0-9.0	55041	33
Klockner		
Micro VEGA 3.0	55074	66
Narrow VEGA 3.5	55065	57
Regular VEGA 4.0, 4.5	55034	26
Narrow KL 3.3	55012	2
Regular KL 3.7, 4.2	55011	1
Wide KL 4.7	55013	3
SK2	55038	30
Essential Cone	55067	59
S3M 3.1	55079	71
Konmet		
Classic 4.0, 4.8	55014	4
Cylindrical 3.5-4.0	55065	57
Cylindrical 4.5-5.5	55057	49
Short Cylindrical 5.0-6.0	55078	70
Conical 4.4/3.4	55065	57
Conical 4.8/3.8, 5.5/4.5	55057	49
Kyocera		
NP FINESIA Bone Level 3.2-3.4	55058	50
RP FINESIA Bone Level 3.7-4.2	55034	26
WP FINESIA Bone Level 4.7-5.2	55014	4
RP FINESIA Tissue Level 3.7-4.7	55014	4
WP FINESIA Tissue Level 4.7-5.2	55014	4
TP/ST POI (HAC) EX 3.7	55085	77
TP/ST POI (HAC) EX 4.2	55025	16
POI (HAC) 2P 3.7/4.2	55004	18
Lasak		
QR BioniQ	55066	58
QN BioniQ	55021	11

	No	Type
Leader		
Internal Hex 3.3-5.5	55034	26
Cono-In	55067	59
Maxillent		
iRaise	55049	41
Medentika		
NI Microcone	55021	11
RI Microcone	55081	73
RI Quattrocone	55081	73
30 AI Quattrocone	55081	73
Procone 3.3	55119	85
Procone 3.8	55031	23
Procone 4.3	55031	23
Procone 5.0	55032	24
Medentis		
BL ICX Premium/Active Master 3.75-4.8	55062	54
Tissue Level ICX Plus 3.45	55039	31
BL ICX Premium/Active Master 3.3	55058	50
Medical Instinct		
Hex 3.0	55071	63
Hex 3.4, 4.0, 5.0	55052	44
Cone 3.0, 3.4, 4.0	55071	63
Cone 5.0	55072	64
Medigma Biomedical		
Mars 3.5-6.0	55057	49
FixTite 3.5-6.0	55057	49
Fix-a-Dent 3.75-6.0	55057	49
Medimecca		
PT SLA Chaorum	55047	39

	No	Type
P-I Brånemark		
HEX 3.3	55115	81
HEX 3.75-5.0	55117	83
AMP 3.3	55039	31
AMP 3.75-4.8	55118	84
MT 3.3-4.8	55046	38
Conical Abutment HEX, AMP MT	55080	72
Paltop		
Standard 3.75-5.0	55034	26
Narrow	55027	19
Conical PCA	55066	58
Park Dental Research		
P series Starvent	55046	38
XP series Starvent	55049	41
Narrow ridge Startanius	55046	38
Phibo		
NP Aurea EVO ☒	55027	19
RP, WP Aurea EVO ☒	55051	43
TSA S3☒	55074	66
TSA S4, S5	55034	26
TSH S3, S4, S5☒	55011	1
PLATON JAPAN		
NC Eight-Lobe Pro 3.5-4.1	55046	38
RC Eight-Lobe Pro 4.5-5.0	55014	4
Type IV PRO	55020	10
Prodent Italia		
Orange Twinner & Prime, Conomet, Collar	55027	19
Fuchsia Twinner & Prime, Conomet, Collar	55051	43
Green Twinner & Prime, Conomet, Collar	55067	59
Yellow Twinner & Prime, Conomet, Collar	55020	10
Blue Prime, Conomet, Collar	55020	10
Qualia		
NeoActive	55046	38
Quest		
NP Hylock	55030	22
RP Hylock	55066	58

	No	Type
Resista		
WIS 3.0-5.0	55049	41
Rex Implants		
TL External hex 1.8	55060	52
TL Internal Hex 2.9	55034	26
BL Internal Hex 2,9	55034	26
Ritter		
Multi Unit	55033	25
Narrow Spiral SB/LA 2.9	55081	73
Standard Spiral SB/LA 3.75	55049	41
ROTT Implants (TRATE)		
R	55065	57
RS Implant		
3.5-6.0	55020	10
Schütz Dental		
Impla cylindrical	55052	44
IMPLA micro Retention cone 4.2	55052	44
SGS Implant		
P1, P7	55049	41
P1D, P5D, P7D	55049	41
P7N	55059	51
SHOFU		
BIOFIX 3.0	55064	56
BIOFIX 4.0	55046	38
BIOFIX 4.5	55046	38
SIC Invent		
MUA	55120	86
SICace 3.4-5.0	55052	44
SICmax, SICtapered 3.7-5.2	55052	44
SICvantage max, SICvantage tapered 3.0	55026	17
SICvantage® max, SICvantage tapered 3.7	55116	82
SICvantage max, SICvantage tapered 4.2-5.2	55054	46
Simply Implants		
Internal Hex	55051	43

	No	Type
Zimmer/Zimvie		
Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41
Screw-Vent, Tapered Screw-Vent 4.7-6.0	55051	43
SwissPlus 3.8	55070	62
SwissPlus 4.8	55014	4
Spline Twist 3.75-5.0	55038	30
NP Eztetic 3.1	55065	57
TSX 3,5	55049	41
TSX 4,5	55051	43
Zuga		
All	55027	19

MULTIPEG assortment

	No	Type
Bego		
Mini Semados	55012	2
S / RI Semados 3.25-5.5	55034	26
RS / RSX Semados 3.0	55082	74
RS / RSX Semados 3.75-5.5	55034	26
SC / SCX Semados 3.25-5.5	55034	26
Bio3		
Advanced and Progressive 3.3-3.8	55025	16
Advanced and Progressive 4.2-5.0	55067	59
BioHorizons		
Grey Laser-Lok 3.0	55034	26
Yellow Internal Laser-Lok 3.5	55049	41
Green Internal Laser-Lok 4.5	55051	43
Blue Internal Laser-Lok 5.7	55051	43
External 4.0	55011	1
Multi-Unit Abutment	55080	72
NP Grey Tapered Pro Conical 3.3-3.8	55066	58
RP Yellow Tapered Pro Conical 4.2-5.2	55066	58
BioMate		
BioMate 3.3	55074	66
BioMate 4.8	55081	73
Biomet/3i/Zimvie		
Purple Certain internal	55036	28
Purple Certain internal 3.4	55036	28
Blue, Yellow, Green Certain internal	55066	58
External Hex	55011	1
Purple T3	55036	28
Blue, Yellow, Green T3	55066	58
Bioner		
HIKELT 3.8-4.7	55034	26
Top DM 3.5-5.0	55046	38
Bionika		
Cortical 4.0	55034	26
CortiLog (CCL) 3.3	55065	57
CortiLog (CCL) 3.8	55023	14
CortiLog (CCL) 5.0	55075	67

	No	Type
Bionnovation		
Biomorse	55023	14
Biotech Dental/ids		
Kontakt, Kontakt-S 3.6-5.4	55045	37
Kontakt 3.0	55130	96
Kontakt-S+ 4.5-5.5	55045	37
Reflect Aspire	55081	73
Reflect Certus	55019	9
Reflect Tapered	55034	26
Reflect Recover	55020	10
Reflect Rapid	55020	10
Biotem		
Yellow M AR 3.7	55066	58
Green R AR 4.2-4.6	55037	29
Blue W AR 5.1	55037	29
Yellow M BR 3.3	55012	2
Green R BR 3.75-4.0	55014	4
Blue W BR 5.0-6.0	55013	3
Yellow M IR 3.5	55028	20
Green R IR 4.1-4.8	55028	20
Blue W IR 4.8-5.5	55028	20
BlueSkyBio		
NP BIO Max All	55046	38
RP BIO Max 4.3-5.0	55037	29
3.7-4.1 BIO Int Hex 3.5	55049	41
4.7 BIO Int Hex 3.5	55051	43
4.7-7.0 BIO Int Hex 4.5	55051	43
Narrow, Regular, Wide BIO One-stage	55014	4
Short, Regular, Wide BIO One-stage	55014	4
Narrow BIO Quattro 3.3	55061	53
Regular BIO Quattro 4.1-5.6	55062	54
3.0 BIO Conus 3.0	55064	56
3.5-4.0 BIO Conus 3.5-4.0	55046	38
BIO Conus 4.5 x 6 mm	55046	38
4.5-5.0 BIO Conus 4.5-5.0	55017	7
BIO Conus 5.0 x 6 mm	55046	38
BIO Trilobe 3.5	55016	6
BIO Trilobe 4.3-6.0	55017	7

	No	Type
GDT		
Regular Internal Hex	55049	41
Narrow Internal Hex	55039	31
Regular Conical connection	55037	29
Narrow Conical connection	55046	38
Geass		
Way Extra	55081	73
Global D		
In-Kone 4.0-5.0	55024	15
WD In-Kone	55133	100
TwinKon 3.5-4.5	55122	88
gmi		
RP avantgard 3.75-4.25	55066	58
WP avantgard 4.75-5.75	55075	67
RP, WP frontier	55051	43
NP phoenix 3.3	55038	30
Std-wide phoenix 3.3-5.0	55011	1
monolith 3.0	55043	35
GP Implants		
SSI Spiral Shape Impant	55051	43
gt medical		
HXE	55011	1
RP HXE	55012	2
WP HXE	55013	3
HXI 4.1	55037	29
OO2 HXI 4.7	55049	41
OCT 4.3	55014	4
Hahn		
Tapered 3.0	55058	50
Tapered 3.5	55046	38
Tapered 4.3	55046	38
Tapered 5.0	55037	29
Tapered 7.0	55048	40
Hi-Tec		
NTI, SLT Tapered Self Thread 3.3, 3.75	55034	26
GTI, WTI Tapered Self Thread	55051	43
TPI Spark	55034	26
Plus Logic	55034	26
LGI Logic 3.0-3.5	55046	38
LGI Logic 4.3-5.0	55037	29

	No	Type
Hiossen		
HS system	55028	20
mini ET system 3.2	55047	39
mini ET system 3.5	55015	5
regular ET system 4.0-7.0	55037	29
single EK system 3.5	55046	38
single EK system 4.0-7.0	55044	36
IDI Evolution		
Standard Implant	55004	18
High B1One	55080	72
Low B1One	55033	25
IDI Implants		
BIO 3.5-5.0	55037	29
ST ID Cam 4.2	55010	13
All ID	55010	13
ID3 3.1	55021	11
Ids/Biotech dental		
Kontakt, Kontakt-S 3.6 - 4.8	55045	37
Aspire Reflectt	55081	73
Certus Reflect	55019	9
Tapered Reflect	55034	26
Recover Reflect	55020	10
Rapid Reflect	55020	10
IML Swiss Dental		
Starfly 3.5-6.0	55034	26
Universe 2.9	55039	31
Universe 3.4-5.0	55024	15
Implance		
Tissue level	55014	4
Regular Bone level	55037	29
Bone level 3.3	55065	57

-

	No	Type
Megagen		
AnyRidge	55051	43
AnyOne Internal	55037	29
Narrow BlueDiamond	55061	53
Regular BlueDiamond	55066	58
Mini	55064	56
MUA	55080	72
Octa Abutment	55121	87
Meisinger Implants		
RP, WP OKTAGON Tissue level	55014	4
NC OKTAGON Bone Level	55061	53
RC OKTAGON Bone Level	55062	54
Microdent		
Universal	55011	1
Microdent 2.8-3.25	55068	60
Microdent 3.3-5.0	55011	1
Microdent 5.5-8.0	55013	3
Genius 3.0	55025	16
Genius 3.5-5.0	55051	43
Ektos 3.2	55023	14
Ektos 3.7-8.0	55051	43
Trylogic 3.5-4.2	55051	43
Trylogic 5.0	55075	67
MIS		
Narrow (C1, V3) Conical Connection	55081	73
Standard (C1, V3) Conical Connection	55087	79
Wide (C1) Conical Connection	55041	33
Narrow Seven and M4	55027	19
Standard Seven and M4	55049	41
Wide Seven 5.0	55051	43
NP, SP, WP Connect Abutment	55115	81
Multi Unit System	55033	25
Mode Implant		
NP Rapid 3,3	55046	38
NP Rapid 3,7	55046	38
RP Rapid 4,1	55027	19
RP Rapid 4,7	55027	19
RP Rapid 5,2	55027	19
NP Shorter 3,7	55046	38
RP Shorter 4,1	55027	19
RP Shorter 4,7	55027	19
RP Shorter 5,2	55027	19
NP Level 3,3	55027	19
NP Level 3,7	55027	19

	No	Type
RP Level 4,1	55027	19
Mode Implant		
RP Level 4,7	55027	19
RP Level 5,2	55027	19
NP Short 3,7	55065	57
RP Short 4,1	55081	73
RP Short 4,7	55081	73
RP Short 5,2	55081	73
myplant GmbH		
myplant two 3.5-4.5	55077	69
myplant bio 4.0	55077	69
Neo Biotech		
IS II 4.0	55020	10
IS III 4.0	55020	10
Neodent		
Morse Connection Cone	55025	16
WS Connection Titamax	55056	48
Connection Grand Morse	55066	58
GM Micro Abutment	55080	72
Helix Short Purple	55056	48
Narrow Connection Grand Morse	55021	11
Zi Ceramic	55019	9
Neoss		
SP 3.5 implants and wider	55004	18
NP 3.25	55027	19
Access Abutment	55060	52

-

	No	Type
SIN		
SW CM Strong	55023	14
SW Int Hex Strong	55035	27
Tryon 3.25-5.0	55011	1
CM Unitite 2.9	55048	76
CM Unitite 3.5	55077	69
CM Compact Unitite 4.0	55085	77
Southern Implants		
Tri-Nex 4.3-5.0	55010	13
Tri-Nex 3.5	55022	12
External Hex	55011	1
MAX Ext Hex	55011	1
0, 12, 24, 36 Co-Axis Ext Hex	55011	1
Piccolo Ext Hex	55050	42
IT-Connection	55014	4
Deep Conical 3.0	55064	56
Deep Conical 3.5-4.0	55046	38
Tapered Deep Conical 5.0	55017	7
Spiral Tech		
Hex connection ESi	55051	43
Hex connection Ultimate	55051	43
Conical ESi 3.5	55044	36
Conical ESi 4.3, 5.0	55037	29
Stark Implants		
Narrow Gravis 3.0-3.5	55064	56
RP Gravis 3.8-6.0	55037	29
Straumann		
Regular & Wide neck Tissue level	55014	4
Narrow Neck Crossfit (NNC) Tissue level	55076	68
Tapered Narrow (NC) Bone level	55061	53
Tapered Regular (RC) Bone level	55062	54
Tapered SC Bone level 2.9	55064	56
BLX all platforms Bone level	55027	19
BLC all platforms Bone Level C	55027	19
TLX all platforms Tissue Level X	55115	81
TLC all platforms Tissue Level C	55115	81
NC, RC Screw retained abutment	55080	72
PURE ceramic 4.1	55047	39

	No	Type
Surgikor		
Fixation 3.0	55058	50
Fixation 3.5	55066	58
Fixation 4.3	55037	29
Solution 3.5, 4.0	55066	58
Solution 4.5, 6.0	55037	29
Immediate all	55049	41
Versatile all	55049	41
Sweden &#038; Martina		
Syra 3.8-5.0	55011	1
Outlink2 3.75-5.0	55011	1
Outlink2 3.3	55079	71
Shelta Power 3.8-5.0	55025	16
Shelta 3.8-5.0	55025	16
Premium POWER 3.3-5.0	55025	16
Premium ONE 3.3-5.0	55025	16
Prama 3.3-5.0	55079	71
CSR-N DAT	55023	14
CSR, DAT connection	55035	27
Shelta Abutment non-EU markets only	55128	94
Youse I	55034	26
Youse V-N	55065	57
Youse V-R	55020	10
Tag Dental		
Axis	55051	43
massif	55051	43
Tatum		
Integrity 3.7-5.0	55015	5
Integrity 6.0-8.0	55040	32
T Implant	55063	55
TBR		
Infinity PF 3.5, 4.0	55004	18
M Implant	55044	36

—

	No	Typ
Brain Base		
B-Type Mytis Arrow 3.7-4.6	55014	4
BS-Type Mytis Arrow 3.7-4.6	55014	4
BW Mytis Arrow 4.7	55014	4
C-Type Mytis Arrow 3.3	55014	4
C-Type Mytis Arrow 3.7-4.6	55020	10
E-Type Mytis Arrow 3.7-4.6	55014	4
EW Mytis Arrow 5.0-7.0	55014	4
Me-Type Mytis Arrow 4.0-5.0	55020	10
St-Type Mytis Arrow 4.0-5.0	55020	10
Bredent		
SKY	55025	16
copaSKY	55027	19
mini2SKY	55033	25
fast & fixed Abutment	55083	75
BTI		
Interna Universal	55051	43
Interna Wide	55051	43
BTK		
EN External Hex	55038	30
ER External Hex	55011	1
EW External Hex	55011	1
IR Internal Hex	55049	41
IM Internal Hex	55051	43
IW Internal Hex	55051	43
KR Hexagonal Morse Taper	55081	73
KW Hexagonal Morse Taper	55056	48
SR Octagonal Morse Taper	55014	4
C-Tech		
EL	55016	6
Camlog		
Screw-Line, Root-Line, Progressive 3.3	55036	28
Screw-Line, Root-Line, Progressive 3.8, 4.3	55031	23
Screw-Line, Root-Line, Progressive 5.0, 6.0	55032	24
iSy	55065	57
Cell		
Cell 3,85	55066	58

	No	Typ
Ceralog		
Hexalobe ceramic	55076	68
Cervico System		
i-Base	55018	8
CLC Scientific		
Conic 3.5-6.0	55081	73
Conelog		
Screw-Line, Progressive 3.3	55066	58
Screw-Line, Progressive 3.8, 4.3	55066	58
Screw-Line, Progressive 5.0	55075	67
Cowellmedi		
Inno Int Octa	55014	4
Inno Sub hex	55037	29
Dentack		
Internal Hex	55056	48
Dental Ratio		
RP, WP Tissue Level	55014	4
NC Bone Level	55061	53
RC Bone Level	55062	54
Dental Tech		
FT3 Implassic 3.25	55039	31
FT3 Implassic 3.75	55049	41
FTP Implassic	55019	9
HX3 Implassic	55011	1
TR2 Implassic	55014	4
TRW2 Implassic	55014	4
Dentalis		
Regular	55051	43

	No	Typ
Implant Direct		
Legacy 1/2/3/4, 3.0/3.2	55012	2
Legacy 1/2/3/4, 3.5/3.7-5.7/7.0	55051	43
RePlant 3.5	55022	12
RePlant 4.3-5.0	55010	13
RePlant 6.0	55032	24
ReActive 3.5	55022	12
ReActive 4.3-5.0	55010	13
3.0 Interactive	55046	38
3.4 Interactive	55037	29
3.0 Simply Iconic	55046	38
3.4 Simply Iconic	55037	29
Screwindirect 5.0 PF 3.7-5.7	55070	62
Implant One		
Std Thread 100 Series	55029	21
Std Thread 200 Series	55029	21
Std Thread (blue) 300 Series	55029	21
Std Thread (purple) 400 Series	55035	27
Std Thread (green) 500 Series	55035	27
Wide Thread (blue) 300 Series	55029	21
Wide Thread (purple) 400 Series	55035	27
Wide Thread (green) 500 Series non-EU markets only	55129	95
Implantswiss		
Tissue level	55014	4
Regular Bone Level	55037	29
Bone Level 3.3	55065	57
Bone Level 2.9	55058	50
Implantvision		
Internal Hex 3.7-5.7	55034	26
Innobiosurg		
FC Magic	55018	8
Magicore 4-6.5	55050	42
Intoss Anchor		
INternal hex 3.7-5.7	55088	80
Intra-lock		
Intrahex 3.4 - 4.5	55051	43
Gold & Blue 3.4 - 6.0	55051	43
InDex 3.4 - 3.75	55015	5
InDex 4.0 - 4.3	55020	10
FlatTop 3.75	55012	2

	No	Typ
IRES		
Narrow External Hex 3.3	55012	2
Regular External Hex 3.7-5.2	55011	1
Narrow Internal Hex 3.3	55087	79
Regular Internal Hex 3.7-5.2	55049	41
Regular Internal Octagon 3.7-5.2	55014	4
Izen Implants		
Narrow I-System X	55064	56
Narrow I-System	55047	39
Mini/Regular/Wide I-System	55020	10
Narrow T-System	55047	39
Mini T-System	55015	5
Regular/Wide T-System	55037	29
Mini/Regular/Wide R-System	55051	43
J Dental Care		
JD Evolution 3.7-6.0	55067	59
Plus+ JD Evolution 3.7-6.0	55067	59
S JD Evolution 3.25	55044	36
S JD ICON ultra 2.75	55064	56
Conical Abutment	55080	72
Kamal Medtech		
i-fix 3,5-5,0	55037	29
s-fix 3,6-6,0	55046	38

	No	Typ
Nobel Biocare		
RP Brånemark System Ext Hex	55011	1
NP Brånemark System Ext Hex	55012	2
WP Brånemark System Ext Hex	55013	3
NP Tri-Lobe Connection	55022	12
RP, WP Tri-Lobe Connection	55010	13
Tri-Lobe Connection 6.0	55032	24
Conical Connection 3.0	55058	50
Narrow Conical Connection	55046	38
Regular Conical Connection	55037	29
Wide Conical Connection 5.5	55048	40
Regular N1	55126	92
Narrow N1	55127	93
NP, RP MUA	55080	72
NP, RP, WP On1	55042	34
Nobel Pearl	55086	78
Noris Medical		
Tuff 3.3-6.0	55034	26
Novem		
3.4	55039	31
3.8	55016	6
4.2, 5.0	55057	49
NucleOss		
Tpure Bone Level	55004	18
T6 Bone Level	55074	66
OCO Biomedical		
Regular Engage 3.25/3.0	55027	19
Regular Engage 4.0/3.5	55051	43
Regular Engage 5.0/4.5	55051	43
Regular Engage 6.0/5.7	55051	43
Engage Plus 4.0-6.0/3.25		
ERI 3.25/3.7	55069	61
ERI 4.0/4.7	55014	4
ERI 5.0/5.7	55014	4
TSI 3.25/3.7	55068	60
TSI 4.0/4.7	55011	1
TSI 5.0/5.7	55011	1
OsseoFuse		
Hexaplus/ HexaPlus-S Sinus 3.75-7.5	55081	73
ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5

	No	Typ
Osstem		
Regular US	55011	1
Mini US	55012	2
Wide US	55013	3
Wide PS US	55011	1
Regular, Wide SS	55028	20
Regular, Regular Short TS 6, 7	55037	29
Mini TS 3.5	55015	5
Mini TS 3.0	55047	39
Osteo Plus		
PF Iridium 3.8, 4.2	55079	71
She 3.45, 3.75	55011	1
Shi 3.75, 4.2	55049	41
Shi 5.0	55051	43
OT Medical		
OT-F2 3.4, 3.8	55053	45
OT-F2 4.1	55073	65
Overmed		
Easyline Easydip	55020	10
ETLINE	55034	26
Oxy Implant		
Internal	55020	10
Narrow Kone Conical	55021	11
Regular Kone Conical	55051	43
Regular PSK Line	55051	43
Mini External	55079	71
Regular External	55011	1
Maxi External	55011	1

	No	Typ
TFI System		
HEX Connection Easy Grip Implants	55087	79
CONE Connection Easy Grip Implants	55067	59
Narrow Connection Easy Grip Implant	55074	66
Thommen		
SPI 3.5	55021	11
SPI 4.0	55065	57
SPI 4.5-6.0	55066	58
ticare		
Mini, Standard Osseus	55011	1
Maxi Osseus	55013	3
Maxi ps Osseus	55011	1
Mini Inhex	55021	11
Standard Inhex	55046	38
Maxi Inhex	55017	7
Standard Osseus Quattro	55011	1
Standard Inhex Quattro	55046	38
Mini Pilar Conico	55021	11
STD Pilar Conico	55046	38
TiTC		
MT IM 3.5 - 5.6	55049	41
TRATE		
All diam ROOTFORM	55065	57
TRI Dental implants		
Tissue level Octa	55014	4
all Tri-Vent	55034	26
Narrow	55025	16
Trinon		
Q2	55011	1
QK	55014	4
QZA 3.5	55049	41
QZA 4.5, 5.7	55051	43
URIS		
Narrow Omni 3.0-3.5	55115	81
Regular Omni 4.0-6.5	55020	10

	No	Typ
Vulkan		
Internal hex	55034	26
NP External hex	55012	2
RP External hex	55011	1
MiP Conical con	55058	50
NP Conical con	55046	38
Winsix		
Internal hex	55020	10
Xgate		
Conical 3.3-5.0	55034	26
Xive		
S-Plus 3.0	55030	22
S-Plus 3.4	55061	53
S-Plus 3.8-5.5	55062	54
TG Plus 3.4-4.5	55018	8
Z-Systems		
Bone level Z5	55131	97
Tissue level Z5	55133	99
Zeramex		
XT Implant 4.2, 5.5	55086	78
Ziacom		
NP Internal hexagonal connection	55027	19
RP Internal hexagonal connection	55057	49
WP Internal hexagonal connection	55057	49
NP External hexagonal connection	55079	71
RP External hexagonal connection	55117	83
WP External hexagonal connection	55014	4
RP Conical connection	55046	38
WP Conical connection	55075	67